

Action Groups guide

English

Action Groups* is a special feature within AirView™ that makes managing your sleep business easy. With Action Groups, patients are sorted into six distinct categories that identify which patients have met compliance guidelines and which ones are struggling with therapy.


To navigate, from the **Patients** menu, select **Action Groups**.

Additional filtering options

Filter groups by clinical user, location, or days since their start of therapy.

Action groups

Patients are sorted into actionable groups based on certain criteria. Click on any of these groups to see patient therapy data for the previous 10 days.



Number of patients

Indicates the number of patients in the group.

Rule description

Hover over any of these groups to view a brief description of the specific rules for each group.

Compliance risk	Ongoing usage follow up	Therapy issues	No data received	Compliance met	Compliance not met
Identifies patients who've demonstrated low device usage prior to achieving initial compliance.	Identifies patients who have demonstrated low device usage after they've achieved initial compliance, or after their first 90 days of therapy.	Identifies patients who've experienced mask leak, AHI or CAI values above the set threshold during a specific time-period.	Identifies patients who have not had their data transmitted for three consecutive days within their first 90 days of therapy, or during the most recent 10-day period after their initial 90 days of therapy.	Identifies patients who have achieved initial compliance based on initial compliance rules.	Identifies patients who have not achieved initial compliance based on initial compliance rules.

*This feature is only available with AirSense™ 10, AirCurve™ 10 and S9 devices.



Action Groups:

Compliance risk		
Rule Based On	Threshold	Timing
Daily Usage Hours in the first 90 days since start of therapy	< 4 hours of usage per night	3 days in the first 7-day period
		5 days in the first 14-day period
		10 days in the first 30-day period
Initial compliance not met at day 45, 60 and 75		

Ongoing usage follow up		
Rule Based On	Threshold	Timing
Daily Usage Hours after the first 90 days since start of therapy	< 4 hours	10 days in a row

Compliance met		
Rule Based On	Threshold	Timing
Initial compliance achieved	> = 4 hours usage per night	21 days out of a 30-day continuous period within the first 90 days since start of therapy
	Based on initial compliance rule of payor template selected	Varies

Therapy issues		
Rule Based On	Threshold	Timing
95th % Mask Leak	> 24 L/min	3 days in a row in the first 90 days
	> 24 L/min	5 days in a row after the first 90 days
AHI	> 10 events per hour	5 days in a row in the first 90 days
	> 10 events per hour	7 days in a row after the first 90 days
CAI	> 5 events per hour	For 3 out of 5 days

No data received		
Rule Based On	Threshold	Timing
No Data Received	N/A	3 days in a row DURING first 90 days
		10 days in a row AFTER initial compliance is achieved or after first 90 days

Compliance not met		
Rule Based On	Threshold	Timing
Initial compliance not achieved	< = 4 hours usage per night	21 days out of a 30-day continuous period within the first 90 days since start of therapy
	Based on initial compliance rule of payor template selected	Varies



ResMed AirView™

Patients Business My profile Administration

Mark as reviewed View reviewed patients

Therapy issues

[Back to Action Groups](#)

Name	Initial compliance	Day	Last 30	01/23	01/24	01/25	01/26	01/27	01/28	01/29	01/30	01/31	02/01
Sample, Patient	✗	9	10%	■	■	■	■	■	■	■	■	■	M
Sample, Patient	✗	9	10%	■	■	■	■	■	■	■	■	■	M
Sample, Patient	✗	9	10%	■	■	■	■	■	■	■	■	■	M
Sample, Patient	✗	9	10%	■	■	■	■	■	■	■	■	■	M
Sample, Patient	✗	9	10%	■	■	■	■	■	■	■	■	■	M
Sample, Patient	✗	9	10%	■	■	■	■	■	■	■	■	■	M

Patient group list

When you select a group, you'll see the list of patients in that specific group. This **Therapy Issues** group lists patients who have mask leak, AHI or CAI values above the set threshold for a specific time period.

Action menu

The Action menu is displayed when you select a specific patient record. It lets you easily access usage charts, view device issues via the **Remote Assist** feature, add notes, or create a report to review patient therapy data.

- With the menu, you can also:
- Access prescription settings
 - Review all logs
 - Mark patient records as reviewed from each group

Name	Initial compliance	Day	Last 30	01/23	01/24	01/25	01/26	01/27	01/28	01/29	01/30	01/31	02/01
Sample, Patient	✗	9	10%	■	■	■	■	■	■	■	■	■	M
Sample, Patient	✗	9	10%	■	■	■	■	■	■	■	■	■	M
Sample, Patient	✗	9	10%	■	■	■	■	■	■	■	■	■	M
Sample, Patient	✗	9	10%	■	■	■	■	■	■	■	■	■	M
Sample, Patient	✗	9	10%	■	■	■	■	■	■	■	■	■	M
Sample, Patient	✗	9	10%	■	■	■	■	■	■	■	■	■	M

In groups

If your patient is in multiple Action Groups, these will be shown here.



View recent logs

From the Action menu, click **View recent logs** to see the rules that placed patients into each Action Group. You can see all Action Group logs when you click **View All Action Group logs**.

Sample, Patient 11 60

02/02/2017 Added to Ongoing usage follow up
Usage < 4 for 10 days in a row

[View all Action Group logs](#)

Ok

View reviewed patients

Click **View reviewed patients** to display all patients that have been marked reviewed in the last 48 hours.

Restore selected to group

Reviewed from Therapy issues

[Back to patients](#)

Name	Initial compliance	Day	Last 30	01/27	01/28	01/29
Sample, Patient	X	9	10%	Green	Green	Green
Sample, Patient	X	9	10%	Green	Green	Green
Sample, Patient	✓	72	53%	Green	Red	Red

Mark as reviewed

Once you've taken action on a patient issue (e.g. called the patient, adjusted comfort settings), you can remove the patient's record from the Action Group. From the **Patients>Action Groups** menu, choose the patient record and click **Mark as reviewed**. You can also select multiple or all patients in a group and mark them as reviewed.

Mark as reviewed > View reviewed patients >

Ongoing usage follow up

[Back to Action Groups](#)

Name
Sample, Patient
<input checked="" type="checkbox"/> Sample, Patient

Restore reviewed patients

Select the patient and click **Restore selected to group** to return a patient to the Action Group.

Restore selected to group

Reviewed from Therapy issues

[Back to patients](#)

3 selected [Select all 3 in a group](#)

<input checked="" type="checkbox"/> Sample, Patient
<input checked="" type="checkbox"/> Sample, Patient
<input checked="" type="checkbox"/> Sample, Patient