



Practical guide ASV – ST – iVAPS in sleep labs

From diagnosis to therapy

Abbreviation key

ABG	arterial blood gas	MIP	maximum inspiratory pressure
ASV	adaptive servo-ventilation	MSA	mixed sleep apnea
BR	backup rate	MV	minute ventilation
CompSA	complex sleep apnea	NIV	non-invasive ventilation
COPD	chronic obstructive pulmonary disease	NH	nocturnal hypoventilation
CSA	central sleep apnea	NMD	neuromuscular disease
CSR	Cheyne-Stokes respiration	OHS	obesity hypoventilation syndrome
EPAP	expiratory positive airway pressure	OSA	obstructive sleep apnea
FEV1	forced expiratory volume	PaCO₂	arterial partial pressure of carbon dioxide
FVC	forced vital capacity	PS	pressure support (IPAP – EPAP)
FIO₂	fraction of inspired oxygen	PSG	polysomnography
iBR	intelligent backup rate	RR	respiratory rate
IPAP	inspiratory positive airway pressure	SaO₂	arterial saturation of oxygen
IV	invasive ventilation	SDB	sleep-disordered breathing
iVAPS™	intelligent Volume-Assured Pressure Support	ST	spontaneous-timed (ventilation mode)
LPM	liters per minute	TECSA	treatment-emergent central sleep apnea

Editorial

Obstructive sleep apnea (OSA) represents the most common form of sleep-disordered breathing (SDB). However, many other more challenging SDB conditions like central sleep apnea (CSA) and various forms of hypo- and hyperventilation exist. CPAP and APAP therapy – gold standard therapies for OSA – may be inadequate in these cases. More sophisticated therapy devices and modes are required.


This practical guide is intended to help you diagnose various SDB conditions and treat them using appropriate therapy modes (e.g. ASV*, ST or ResMed iVAPS†) and devices in accordance with existing medical guidelines. It is built on a simplified, clearly laid out presentation that takes into account both medical practice and proven techniques. For more information, including the US Centers for Medicare and Medicaid Services guidelines, visit CMS.gov.¹

* ASV therapy is contraindicated in patients with chronic, symptomatic heart failure (NYHA II–IV) with reduced left ventricular ejection fraction (LVEF \leq 45%) and moderate to severe predominant central sleep apnea

† iVAPS therapy mode is indicated for patients weighing 30 kg (66 lbs) and above

¹ US Centers for Medicare & Medicaid Services. CMS.gov: License Agreements. Baltimore, MD. Accessed online on March 20, 2019: www.cms.gov/medicare-coverage-database/details/lcd-details.aspx?LCDId=33800&ver=13&CoverageSelection=Local&ArticleType=All&PolicyType=Final&ts=All&CptHcpcsCode=e0601&bc=gAAAAACAAAAA&

Indications and therapy modes




Hypocapnia versus hypercapnia

	Hypocapnia	Hypercapnia
Respiration	= Hyperventilation (breathing is too fast and/or too deep)	= Hypoventilation (breathing is too slow and/or too shallow)
Values	<ul style="list-style-type: none"> • $\text{PaCO}_2 < 35 \text{ mmHg}^*$ 	<ul style="list-style-type: none"> • $\text{PaCO}_2 \geq 45 \text{ mmHg}^\dagger$ • $\text{SaO}_2 \leq 88\%$ for \geq a cumulative 5 minutes, minimum 2 hours • $\text{FEV1/FVC} \geq 70\%$
Pathophysiological cause	<ul style="list-style-type: none"> • Hypersensitivity of the respiratory center with minor fluctuations in PaCO_2, therefore normocapnia to hypocapnia is stated by blood gas analysis • PaCO_2 below apnea threshold during sleep > reactive central apnea • PaCO_2 above apnea threshold during sleep > excessive restarting of respiration = hyperventilation 	<p>Reduced gas exchange through:</p> <ul style="list-style-type: none"> • Hypoventilation (in association low tidal volume, inadequate ventilation) • Ventilatory insufficiency (in association with OHS, COPD, thoracic-restrictive or neuromuscular diseases)
Effect	Increased sensitivity of the respiratory center to changes of PaCO_2	High PaCO_2 /low SaO_2
Therapy target	<ul style="list-style-type: none"> • Stabilize breathing • Reduce pronounced PaCO_2 fluctuations and CSA stabilizes periodic breathing 	<ul style="list-style-type: none"> • Improve gas exchange (by supporting ventilation)
Therapy result	<ul style="list-style-type: none"> • Stabilized PaCO_2 • Ideally higher PaCO_2 level towards normocapnic range 	<ul style="list-style-type: none"> • PaCO_2 reduction • Increased SaO_2

[†]Respiratory Assist Device (RAD) Qualifying Guidelines. CMS.gov.

^{*}Refer to full reimbursement guidelines on page 7

From diagnosis to therapy



*Per RAD Qualifying Guidelines. CMS.gov. For CSA diagnosis, central apnea–central hypopnea index (CAHI) is defined as the average number of episodes of central apnea and central hypopnea per hour of sleep without the use of a PAP device.

†ASV therapy is contraindicated in patients with chronic, symptomatic heart failure (NYHA II–IV) with reduced left ventricular ejection fraction (LVEF $\leq 45\%$) and moderate to severe predominant central sleep apnea

Hypercapnia

Diagnosis & indication*

- **COPD:**[†] After initial treatment with E0470 device, ABG shows PaCO₂ worsens ≥7 mmHg compared to original ABG and PSG demonstrates SaO₂ ≤88% for ≥ a cumulative 5 minutes, minimum 2 hours nocturnal recording time while on an E0470 and AHI < 5 OR no sooner than 61 days after initial use of E0470, ABG shows PaCO₂ ≥ 52 mmHg and sleep oximetry demonstrates SaO₂ ≤88% for ≥ a cumulative 5 minutes, minimum 2 hours nocturnal recording time.
- **Restrictive:**[†] Perform one of the following: Either ABG with PaCO₂ ≥ 45 mmHg OR SaO₂ ≤ 88% for ≥ 5 minutes for at least 2 hours nocturnal recording time OR FVC < 50% of predicted or MIP < 60 cmH₂O (NMD only).
- **Hypoventilation:**[†] (After initial treatment with standard bilevel) FEV1/FVC ≥ 70% and PaCO₂ worsens ≥ 7mmHg compared to initial ABG OR A facility-based PSG or HST demonstrates SaO₂ ≤ 88% for ≥ 5 minutes of nocturnal recording time (minimum recording time of 2 hours) that is not caused by obstructive upper airway events – i.e. AHI less than 5 while using an E0470 device.

- NMD (reduced or shallow ventilation in sleep)
- COPD and OSA (overlap syndrome)
- Additional indications for ventilatory support

Therapy approach & target

First and always

- **Provide PS (IPAP – EPAP)** to reduce PaCO₂
- **Backup rate** to ensure adequate minute ventilation and to relieve work of breathing if necessary

Heavily fluctuating PaCO₂ or OSA dependent on position or sleep stage

- **Provide required volume with sufficient PS** to reduce PaCO₂
- **iBR** for better patient-device synchrony

Therapy devices & modes

ST

iVAPS[‡]
iBR

Advanced
Ventilation



AirCurve™ 10

ST-A



Stellar™ series



Astral™ series






Not suitable for periodic breathing (or alternating hypoventilation/hyperventilation)

* The reimbursement information is being provided on an "as is" basis with no express or implied warranty of any kind and should be used solely for your internal informational purposes only. The information does not constitute professional or legal advice on reimbursement and should be used at your sole liability and discretion. All coding, coverage policies and reimbursement information are subject to change without notice. ResMed does not represent or warrant that any of the information being provided is true or correct and you agree to hold ResMed harmless in the event of any loss, damage, liabilities or claims arising from the use of the reimbursement information provided to you. Before filing any claims, it is the provider's sole responsibility to verify current requirements and policies with the payer.

[†] Per RAD Qualifying Guidelines. CMS.gov.


[‡] iVAPS therapy mode is indicated for patients weighing 30 kg (66 lbs) and above

ResMed therapy device modes

	ASV / ASVAuto	ST	iVAPS* iBR	
Functions	<div></div> <ul style="list-style-type: none">• PS: automatically provides 90% of MV (last 3 minutes)• RR: automatic BR• EPAP: manual in ASV mode• AutoEPAP: set min/max range in ASVAuto mode	<div></div> <ul style="list-style-type: none">• IPAP: manual• EPAP: manual• RR: manual	<div></div> <ul style="list-style-type: none">• PS: automatic for patient's alveolar ventilation target• RR: automatic BR with target patient rate (iBR)	
Information for prescription	<ul style="list-style-type: none">• ASV confirmed in titration• PS (min./max.)• EPAP or AutoEPAP (min./max.)	<ul style="list-style-type: none">• IPAP/EPAP• BR• TiControl™: Ti Min, Ti Max, Ti• Trigger/Cycle	<ul style="list-style-type: none">• iVAPS*• PS (min./max.)• EPAP• Target patient rate for iBR• Ti Min/Ti Max	
Therapy device	<div>AirCurve 10 ASV<ul style="list-style-type: none">• Up to 25 cm H₂O</div>	<div>AirCurve 10 ST<ul style="list-style-type: none">• Up to 25 cm H₂O</div>	<div>AirCurve 10 ST-A<ul style="list-style-type: none">• Up to 30 cm H₂O</div>	<div>Stellar series<ul style="list-style-type: none">• Up to 40 cm H₂O</div>

* iVAPS therapy mode is indicated for patients weighing 30 kg (66 lbs) and above

Practical tips



Minimize
leak

Mask adjustment or change, chin strap, full face mask

- Assures mode functionality
- Increases comfort and compliance
- Add/adjust humidification



Observe
therapy

Regular re-evaluations by ResMed AirView™* or ResScan™

- Confirmation of therapeutic efficacy and compliance
- Adjust TiControls and synchrony features

ResMed offers tailored treatment options to sleep labs.
For more information about our products, visit [ResMed.com](https://www.resmed.com).

* When a ResMed device is paired with the ResMed Connectivity Module (RCM), wireless connectivity is enabled, allowing key patient therapy data to be transmitted directly from the device to the ResMed secure, cloud-based management system, AirView

Notes

Notes



Distributed by ResMed Corp, 9001 Spectrum Center Boulevard, San Diego, CA 92123 USA +1 858 836 5000 or 1 800 424 0737 (toll free). See [ResMed.com](https://www.resmed.com) for other ResMed locations worldwide. AirCurve, AirView, Astral, iVAPS, Stellar, ReScan and TiControl are trademarks and/or registered trademarks of the ResMed family of companies. Specifications may change without notice. For patent and other intellectual property information, see [ResMed.com/ip](https://www.resmed.com/ip). © 2019 ResMed. 10111169/1 2019-06