



## Review Patients From Groups

### **Description:**

After actioning a patient within an Action Group, you can mark the patient as reviewed to remove them from the group.

### **To review patients from an Action Group:**

1. While in an Action Group, select the check box next to a patient's name.  
*Note: You can review multiple patients from the group at the same time by selecting more than one check box.*
2. From the toolbar, click **Mark as reviewed**.
3. If applicable, select additional groups to review patient from. Click **Mark as reviewed** again to confirm.

### **To review one patient from an Action Group while adding a note:**

1. While in an Action group, click a patient's name.
2. Click **Add notes** and enter relevant patient notes.
3. Click **Save and mark as reviewed**.
4. The patient will be reviewed from the group and the note will be saved to their AirView profile.

### **To restore reviewed patients:**

Each patient stays in the restored patients list for 48 hours.

1. From the patient group screen, click **View reviewed patients**.
2. Select the check box next to one or more patient's names.
3. Click **Restore selected to group**.