

Deleting a Modem

Quick Reference Guide



Why?

Perform this procedure if you want to:

- Delete the modem from one patient and assign it to another.
- Replace the patient's current modem.

Note:

If you use any healthcare applications other than AirView™, you must also make the changes in that system.

How?

- 1 Sign in to U-Sleep.
- 2 Click the **People** tab.
- 3 Search for the patient that you want to delete the modem from.
- 4 Click the **Monitoring** tab.
- 5 Click **Edit** to place the tab in edit mode.
- 6 In the **CPAP Monitoring** section, click **Edit** to access the patient's current monitoring details.
- 7 Highlight and delete the modem serial number. (To replace the patient's modem, enter the new serial number.)
- 8 Click **Save**.

U-Sleep - Powered by Umbian Compliance Engine

Edit Save Delete Close Help

Device Manufacturer: * ResMed

Transmission: * AirView

Model: * S9 Elite

Organization: Sleep Inc

Location: San Diego

Clinical User: asmith

CPAP Serial Number: * 34590327523

Modem Serial Number: 12345678910

Data required by: * 11 : 00 PM

Setup date: (mm-dd-yyyy) * 10-21-2014

Begin notification: (mm-dd-yyyy) * 10-21-2014

* Patient identity and data are fictional.

Distributed by ResMed Corp, 38 Solutions Dr., Halifax, NS B3S 0H1 Canada. See www.resmed.com for other ResMed locations worldwide. U-Sleep, AirView and S9 are trademarks and/or registered trademarks of the ResMed family of companies. © 2017 ResMed. RH-102014/2 2017-03

U-Sleep Support

support@u-sleep.com | 1-800-424-0737

