

Using the Action View Dashboard

Quick Reference Guide



Action View

To use Action View, you must ask a U-Sleep™ representative to set you up. After you are set up, Action View displays by default when you sign in. This view automatically groups patients so you can easily determine what action needs to be taken, streamlining your workflow.

Time-based groups
View patients based on their length of time on therapy.

Usage-issue groups
View patients who are at risk, based on their device usage.

Therapy-issue groups
View patients who are struggling, based on key therapy metrics.

Compliance groups
View patients who have met compliance and those who have not.

Three-bar menu
Click this menu to:

- Change locations
- Switch views
- Filter patients
- Add patients

Long-term monitoring group
View patients who are at risk based on key therapy metrics after their initial 90 days on therapy.

Time-based groups	Usage-issue groups	Therapy-issue groups	Compliance groups
1 patient 7 Day All Patients	1 patient 7 Day At Risk	9 patients No Data Transmitted	1 patient Payor Compliance
0 patients 30 Day All Patients	1 patient 14 Day At Risk	5 patients Therapy Issues	0 patients Payor Compliance Not Met
1 patient 60 Day All Patients	0 patients 30 Day At Risk	1 patient High CAI	
0 patients 90 Day All Patients	1 patient 45 Day At Risk	0 patients Post 90 Day Monitoring	

© 2010-2014 Umbian Inc. | Terms of Service | Privacy Policy | Umbian.com | U-Sleep.com | Powered by Umbian Compliance Engine

Using the Action View Dashboard

Quick Reference Guide



Action view details

If you click a group on the **Action View** page, U-Sleep shows you the details.

View patients

See a list of the patients who require your attention based on the group you clicked.

Common tasks

Click a patient's name to see some basic information and a listing of common tasks.

Notes

Click a Notes icon to see the note. U-Sleep shows these icons for any notes that you added for a patient within the last 10 days.

The screenshot shows the U-Sleep dashboard interface. At the top, it says "Signed in as: Emma Smith" with links for "My Profile", "Help", and "Sign Out". Below the navigation bar (HOME, LOCATION, PEOPLE, REPORTS), the main heading is "Therapy Issues". There are filters for "Patient reviewed" and "View reviewed patients", along with "Export to Excel" and "Export to PDF" buttons. A search box labeled "Enter last name" and a "Find" button are also present. The main content is a table with columns for "Patient", "Compliance", "Day", and a date range from "24 Oct" to "3 Nov". The table lists four patients: Arbuckle, Laura (Compliance: ✓, Day: 63), Abdulla, Mohamed (Compliance: --, Day: 55), Li, Wei (Compliance: ✓, Day: 21), and Zhang, Jing (Compliance: --, Day: 21). Each row has a grid of colored squares representing daily compliance status and a checkbox on the right. A blue pencil icon is visible next to the "Li, Wei" row.

Payor compliance status

These icons show you if the patient met initial payor compliance. To learn more about these icons, see the next page.

Day

See the number of days since initial device setup.

Daily compliance status

These icons show you the patient's compliance status for each day. To learn more about these icons, see the next page.

Using the Action View Dashboard

Quick Reference Guide



Compliance status icons

Payor compliance status	
✓	Compliant
✗	Non-compliant
--	Compliance undetermined

Top section indicates
AHI status

Middle section indicates
overall compliance

Bottom section indicates
leak status



Daily compliance status			
	Compliant		Non-compliant
	Compliant—high AHI		Non-compliant—high AHI
	Compliant—high leak		Non-compliant—high leak
	Compliant—high AHI, high leak		Non-compliant—high AHI, high leak
	Compliant—no AHI data		Non-compliant—no AHI data
	Compliant—no leak data		Non-compliant—no leak data
	Compliant—high AHI, no leak data		Non-compliant—high AHI, no leak data
	Compliant—no AHI data, high leak		Non-compliant—no AHI data, high leak
	Compliant—no AHI data, no leak data		Non-compliant—no AHI data, no leak data
	Not monitoring—no data available		Device not used

Distributed by ResMed Halifax, ResMed 38 Solutions Dr., Halifax, NS B3S 0H1 Canada +1 (877) 242 1703 Ext 369235. See www.resmed.com for other ResMed locations worldwide. U-Sleep is a trademark and/or registered trademark of the ResMed family of companies. © 2016 ResMed. RH-102026/1 2016-04

U-Sleep Support

support@u-sleep.com | 1-800-424-0737

