



## View and Filter Home Sleep Test Patients Page

### Description:

The Diagnostic patients screen displays all home sleep test patients within your organization.






### To view patients list:

1. From the **Patients** menu, select **Diagnostic**.
2. Select the patient name from the list.

**Note:** The status icon changes to red if a patient has not moved to the next status after 7 days. When the next step is completed, the icon changes back to green.

### To filter patients list:

1. From the **Patients** menu, select **Diagnostic**.
2. Click **Show Only**.
3. Select the required criteria from the drop-down lists.
4. Click **Apply**.
5. Each group will be filtered to only display patient records that match the selected criteria.
6. To remove the filter, unselect the checkbox.

Status	Description
	<b>No device</b> No icon is shown when a patient has just been set up with no device association.
	<b>Patient testing</b> When a device has been associated and awaiting a home sleep test.
	<b>Physician reviewing</b> Patient data has arrived and is awaiting review and interpretation from the physician.
	<b>Results complete</b> Physician has reviewed patient data and interpretation has been saved.
	<b>Awaiting prescription<sup>^</sup></b> Interpretation has been completed and awaiting prescription for therapy.
	Prescription completed and can be on several stages ( <b>Ready to send, Results sent, Patient accepted or Patient returned</b> ). <sup>^</sup>

Feature not available in all regions