



Update Physician Access

Description:

All Physicians associated with an organization can view patient data, add notes and (depending on their permission level), can edit therapy settings. You can adjust the security settings that apply to all associated physicians.

To edit an organization's security details:

1. From the **Administration** menu, click **Organization** details.
2. Select **Security** from the left side of the screen.
3. Click **Edit** and select the preferred settings:
 - *View patient data, add notes only*
 - *View patient data, add notes and edit device settings.*
4. Click **Save**.

Note: The security settings applied to an organization apply to all physicians associated with the organization to access patient data